



ZBJ01 Intelligent Disintegration Apparatus

INTELLIGENT DISINTEGRATOR



Security



Intelligence



宁波新芝生物科技股份有限公司
NINGBO SCIENTZ BIOTECHNOLOGY CO., LTD

Product Description

- With a 3-cup design, 3 sets of disintegration time limit tests can be performed simultaneously.
- The overall water bath design is adopted, and the temperature of the water bath is more stable and uniform.
- Using high-precision platinum resistance Pt1000 sensor, the temperature is more accurate.
- The use of a quick-disassemble gondola card can quickly disassemble and assemble the gondola.
- Color touch screen operation, automatic positioning of sandblasting bottom plate.
- 3 cups independent lighting system for easy observation.
- It has the functions of audit trail and authority management.

Features

Color touch screen

DC brushless circulating pump

One-piece paint injection shell

Temperature protection switch

With positioning sandblasting base plate

Audit trail

Integrated water bath

Authority management

Quick disassembly and assembly of gondola cards

High-precision platinum resistance Pt1000 sensor

Independent lighting system

Technical Parameter

Number of gondolas	3
Lifting distance	55mm±1mm
Upgrade frequency	31 times/min
water bath temperature	Room temperature~45°C
temperature resolution	0.01°C
Temperature control accuracy	37°C±0.3°C
Screen aperture	Standard 2mm; optional 1mm, 0.71mm, 0.42mm
authority management	Three levels, seven users
Audit trail	96 pieces ×366 days
Login Method	Username + Password, Fingerprint
print online	Support (optional)

